## Personal Services - High Level Disinfectant Log

Reusable high level disinfectant (HLD) must be tested daily to confirm its minimum effective concentration (MEC). This log sheet can help you to meet this requirement. Always follow the disinfectant and test strip manufacturer's instructions. Use one log sheet for each reuse period.

## **High Level Disinfectant Info**

Name of Disinfectant	
Lot #	
Expiry Date	
Date Opened / Activated	
Date Poured	
# of Reuse Days	
Prepared By	

## **HLD Test Strip Info**

Name of Test Strip				
Lot #				
Expiry Date				
Date Opened				
Test Result	Pass <b>√</b>	Fail X	Pass <b>√</b>	Fail X
Tested By				
	*Use the 2nd column i	if you open a new	bottle of test strips d	uring the reuse period

## **Minimum Effective Concentration Test Log**

Day	Date		Test Result				 ]	T	Test Result			$\top$		T	Test Result		
		Time	Pass V	Fail X	Tested By	Day	Date	Time	Pass V	Fail X	Tested by	Day	Date	Time	Pass V	Fail X	Tested by
1						11						21					
2						12						22					
3						13						23					
4						14						24					
5						15						25					
6						16						26					
7						17						27					
8						18						28					
9						19						29					
10						20						30					

©2020 Alberta Health Services, Safe Healthy Environments, ahs.ca/eph



